


Large Incubators

Posted By [PePressiCon](#) On 25. November 2008 @ 09:06 In | [Comments Disabled](#)

Description

Large Chamber Incubation is the best solution for complex long-term experiments because it guarantees a constant temperature of all devices inside the incubator, including large parts of the microscope. Possible temperatures range from ambient temperature up to 40°C. A CO₂ and O₂ Module may also be added as an upgrade, if required. A variety of components can be used together with one of the large incubators.

Schematic overview



Article printed from PeCon GmbH | Microscope Incubators:
<http://www.pecon.biz>

URL to article: http://www.pecon.biz/?page_id=776

Copyright © 2008 PeCon GmbH. All rights reserved.